


FSIS STEC Guidance:

**Testing solutions to meet evolving
regulations and consumer demand**

Michael Eastwood, Ph.D., BioMonitoring Technology Specialist
Justyce Jedlicka, Food & Beverage Regulatory Liaison and Field
Marketing
May 11, 2023



The life science business of Merck KGaA,
Darmstadt, Germany operates as
MilliporeSigma in the U.S. and Canada.

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada.



Disclaimer

- The information in this seminar is given for the purposes of education and discussion.
- It is not intended to be, and it should not be used as a substitute for regulations or regulatory guidance.
- Decisions and actions should be based on the relevant regulations, guidance document and USDA and FDA chapters.
- Statements and opinions expressed are of the presenter and are not necessarily the views of MilliporeSigma.



01

regulatory Landscape and BACKGROUND

Justyce Jedlicka, Food & Beverage Regulatory Liaison and Field
Marketing

Millipore®

Preparation, Separation,
Filtration & Monitoring Products

Background

February 2023

FSIS implementation of new STEC testing

November 2022

FSIS announcement of cessation of N60 excision sampling

January 2023

FSIS announced expansion to additional raw beef products

November 2022

FSIS announcement of STEC testing and open comment period



STEC - Shiga Toxin Producing E. coli

STEC Facts 2015 - 2019



12,866

Infections

- 22% hospitalized
- 15% travel-associated
- 5% outbreak-associated
- <1% died



July

Highest Month of incidents

- July trends 13.3% cases
- >10% infection rate occurs June - August



O157

Largest identified serotype

- Top 5 serotypes include
- CIDT+(not serogrouped) - 41.43%
- O157 - 17.44%
- O26 - 8.34%
- O103 - 8.29%
- Non-O157, not serogrouped - 6.96%



Minnesota

State with largest rate of infection

- Top 8 States
- Minnesota - 9.0
- California - 7.68
- Oregon - 6.78
- Colorado - 6.7
- Tennessee - 4.66
- New Mexico - 4.62
- Connecticut - 4.11
- Georgia - 3.66

FSIS Directive 10010.10 Revision 5



FSIS Directive 10,010.1 Dashboard | *Sampling Verification Activities for Shiga Toxin-Producing Escherichia Coli (STEC) in Raw Beef Products*

TASK CODE

MT05 - Raw Ground Beef or Veal Sampling - Retail

[Search posted Q/As with "askFSIS" in the title](#)

[Submit question to askFSIS](#)

PREPARE	<p>SCHEDULE THE SAMPLING TASK:</p> <ul style="list-style-type: none"> Select correct assignment Directive 13,000.2 Filter for "Lab Sampling" Consider time needed for establishment to hold lot Find sample task on list, and select "Add" Schedule collection date 	<p>RECORDS TO COMPLETE:</p> <ul style="list-style-type: none"> * Notify plant management before collecting sample (i.e., hold sampled lot) Lab Form Questionnaire as part of PHIS Task <ul style="list-style-type: none"> Data specific to sample collection Supplier information (in case sample is positive) Directive 10,010.3 Gather source material information when collecting a sample.
PERFORM	<p style="text-align: center;">MT05 - Raw Ground Beef or Veal Sampling - Retail</p> <p style="text-align: center; font-size: 0.8em;"><i>Each MTxx Sampling Task is based on the "source material" and "finished product"</i></p> <div style="display: flex; justify-content: space-between; align-items: flex-start;"> <div style="width: 30%; text-align: center;"> <p>SOURCE MATERIAL</p>  </div> <div style="width: 35%; text-align: center;"> <p>PROJECT SUMMARY</p> <p>Retail Firm (Compliance)</p> <hr/> <p>Ground beef from any source, sampled in the final package.</p> </div> <div style="width: 30%; text-align: center;"> <p>SAMPLE METHOD</p>  <p style="font-size: 1.5em; font-weight: bold;">2-lb Total</p> </div> </div> <p style="text-align: right; font-size: 0.8em;">N60 Sampling Video Cloth sampling video</p>	
RESPOND	<p>COMPLETE THE SAMPLING TASK:</p> <ul style="list-style-type: none"> Complete Lab Form Questionnaire Verify address in PHIS and on shipping label match Complete PHIS Task File supplier information for sample Verify product is held pending adulterant test result Wait for sample results (typically 2-3 days) 	<p>ACTIONS IN RESPONSE TO STEC PRESUMPTIVE POSITIVE:</p> <ul style="list-style-type: none"> Review supplier and lotting information Verify all affected product is on hold until result is "confirmed" Await confirmation <ul style="list-style-type: none"> If confirmed positive, verify disposition If confirmed positive, verify HACCP requirements Directive 10,010.2



USDA vs EU STEC Testing

US

- O26
 - O45
 - O103
 - O111
 - O121
 - O145
 - O157
-
- FSIS Directive 10010.1 - Sampling Verification Activities for Shiga Toxin-Producing Escherichia Coli (STEC) in Raw Beef Products - Revision 5
 - BAM Chapter 4A: Diarrheagenic Escherichia coli
-
- Upcoming:
 - FDA Leafy Greens Action Plan
 - USDA 2025 Roadmap




EU

- O26
 - O103
 - O111
 - O145
 - O157
-
- ISO/TS 13136 - Microbiology of food and animal feed — Real-time polymerase chain reaction (PCR)-based method for the detection of food-borne pathogens — Horizontal method for the detection of Shiga toxin-producing Escherichia coli (STEC) and the determination of O157, O111, O26, O103 and O145 serogroups
-
- Upcoming:
 - New Method
 - Detection of STEC
 - Serogroup / Virulence factor



FDA Leafy Greens Action Plan



The graphic features a blue background with a white wavy line at the top. The FDA logo and 'U.S. FOOD & DRUG ADMINISTRATION' are in the top right. The title 'Leafy Greens STEC Action Plan' is in large white font. Below are three columns of numbered items under the headings 'Prevention', 'Response', and 'Addressing Knowledge Gaps'. A photo of a man with a mustache, wearing a plaid shirt and a dark vest, kneeling in a field and holding a head of lettuce, is in the bottom right. At the bottom left, there is a hand icon pointing to a link and two URLs. The date 'April 2021' is in the bottom right.

FDA U.S. FOOD & DRUG ADMINISTRATION

Leafy Greens STEC Action Plan

Prevention	Response	Addressing Knowledge Gaps
1. Advance Agricultural Water Safety	8. Investigation Reports	15. Longitudinal Studies
2. Enhance Inspections, Audits and Certification Programs	9. Conduct Follow-Up Surveillance During the Fall 2020 California Growing/ Harvest Season	16. Data Mining and Analytics on Previous Outbreaks
3. Buyer Specifications	10. Promote Tech-Enabled Traceability	17. Adjacent and Nearby Land Use
4. Leafy Greens Data Trust	11. Improve Utilization of Shopper Card Data	18. Compost Sampling Assignment with California
5. Microbiological Surveys for STEC Detection and Enhanced Sampling Protocols	12. Accelerate Whole Genome Sequencing Data Submissions by States	
6. Increase Awareness and Address Concerns Around Adjacent and Nearby Land	13. Advance Root Cause Analysis Activities	
7. Establish and Strengthen Regular Outreach and Communication Programs for Stakeholders in Growing Regions	14. Enhance Outbreak and Recall Communications	

For more information about the Leafy Greens STEC Action Plan, click any link above or visit:
<https://www.fda.gov/food/foodborne-pathogens/leafy-greens-stec-action-plan>
<https://www.fda.gov/food/foodborne-pathogens/leafy-greens-stec-action-plan-accomplishments>

April 2021



USDA Science Blueprint – 2025 Roadmap

Program Themes

1. Sustainable Agricultural Intensification
2. Agricultural Climate Adaption
3. Food & Nutrition Translation
4. Value-Added Innovations
5. Agricultural Science Policy Leadership


Key Programs / Technology Development

1. Poultry: Salmonella & Campylobacter
2. Beef: STEC & Salmonella
3. Deli Meat: Listeria
4. Pork: Salmonella & Yersinia
5. Hemp
6. Expediting testing confirmation methods: WGS & MALDI



Organization

Hierarchy with five segments



Measure SP-1.1.2.1 Reduction in the proportion of poultry samples with Salmonella serotypes commonly associated with human illness

FSIS is analyzing comments for final rule stating intentions to expand non-O157 Shiga Toxin-producing Escherichia coli (STEC) testing to all beef products that are currently sampled for O157:H7

Agency continues to prioritize ready-to-eat (RTE) sampling for Salmonella and Listeria monocytogenes (Lm)

FSIS to evaluate whole genome sequencing complementary technology as part of initiative to Adopt Innovative Approaches and Technologies

FSIS will conduct further in-depth assessments of the sampling and testing protocols identified through the Strategic Assessment of Sampling Resources



02 Effective pathogen testing solutions: ASSURANCE[®] GDS STEC

Michael Eastwood, Technology Specialist

Millipore[®]

Preparation, Separation,
Filtration & Monitoring Products

Accurately Screening For Top 7 STEC is More Difficult than *E. coli* O157

E. coli
O157

PCR designed
for O157-
specific targets

Most positives
can be
culturally
confirmed

Top 7
STEC

Cannot achieve
specificity with
PCR design:
cross-reactors
inevitable

Minority of
positives are
culturally
confirmed*


* Bosilevac & Koohmarie 2012, App. Env. Micro



USDA Criteria for Top 7 STEC



Immunological




Genetic




Conventional PCR Methods Do Not Detect if *eae* and *stx1/2* are in the Same Cell

Step 1. Enrich potential pathogens



Conventional PCR Methods Do Not Detect if *eae* and *stx1/2* are in the Same Cell

Step 2. Eyschepotentialdipl...
Eyschepotentialdipl...



Step 3. Load PCR Tube




Positive




Conventional PCR Methods Do Not Detect if *eae* and *stx1/2* are in the Same Cell

Step 2. Eyschelpotentialipgkphs



Step 3. Load PCR Tube



False Positive



Assurance® GDS

What Is GDS?




Innovative PCR-based molecular pathogen screening method for food-borne pathogens

Developed to meet the daily rigors & requirements specific to food testing labs – adjustable throughput for small or large testing volumes

Designed to provide the most accurate and reliable results – all kits are AOAC validated across a wide range of matrices




Assurance[®] GDS Workflow Fundamentals





Assurance[®] GDS How are we different?

PickPen[®] IMS




Proprietary “intrasolution” **Immuno Magnetic Separation (IMS)** using the PickPen[®]





Assurance[®] GDS DNA Target Amplification and Detection

DNA Amplification
with Lyophilized
PCR Mix



Specific Detection via
fluorescent probe



Automatic Analysis Software



Integrated Internal Control (prevent FN)



Results

- Positive
- Negative
- No Amp



Assurance® GDS STEC Suite

Screening and Serogroup Identification Tools

Step 1: Top 7 STEC Screen

- PickPen® IMS selection of Top 7 STEC
- Differentiation of Big 6 STEC and O157

Step 2: Molecular Serogroup Identification

- Specificity through **both** serogroup and *eae*-subtype PCR targets



Assurance® GDS STEC Suite Screening and Serogroup Identification Tools

Enrich
mEHEC®
Medium
10 Hours



Assurance® GDS STEC Suite


MPX Top 7 STEC Kit Screens for Virulence Genes in Top 7 STEC



Enrich
375g beef
in
mEHEC®
for 10h




IMS of
Top 7
STEC



Specific
Detection
of
Virulence
Genes


Immunological

Genetic



Assurance[®] GDS STEC Suite Conventional PCR Methods Do Not Detect if *eae* and *stx1/2* are in the Same Cell

Step 2. Eyschelpotential plasmids




Step 3. Load PCR Tube




False Positive




Assurance[®] GDS STEC Suite PickPen[®] IMS Enriches for Top 7 STEC




Step 1. Capture pathogens with PickPen[®]




Assurance[®] GDS STEC Suite PickPen[®] IMS Enriches for Top 7 STEC






Step 1. Capture pathogens with PickPen[®] IMS
Step 2. Separate PCR target pathogens



Assurance® GDS STEC Suite Screening and Serogroup Identification Tools



Assurance® GDS STEC Suite MPX ID STEC Identification Kit



Assurance[®] GDS STEC Suite

MPX ID STEC Identification Kit




Big 6
Positive
Sample



IMS of
Big 6
STEC



- *eae1* or ϵ
- *eae2* or θ
- *eae3* or β
- O103
- O121
- O45
- O26
- O111
- O145



STEC ID





Name	Top 6 STEC	O26	O45	O103	O111	O121	O145
O103, epsilon	Positive			+			
O121, epsilon	Positive					+	
O45, epsilon	Positive		+				
O103, beta	Positive			+			
O145, beta	Positive				+		+
O111, theta	Positive				+		
O26, beta	Positive	+					
NC	Negative						

ID of O group
simplifies isolation



Assurance[®] GDS STEC Suite EHEC ID *E. coli* O157:H7 Identification Kit



Assurance[®] GDS STEC Suite **Summary**

