


The new way of making media

The ReadyStream[®] system,
a new solution for culture media
preparation and dispensing




The Life Science business of Merck KGaA,
Darmstadt, Germany operates as
MilliporeSigma in the U.S. and Canada.

Millipore[®]

Preparation, Separation,
Filtration & Monitoring Products

How to get sterile media

Your current options



Standard
Autoclave

Media
Prep

Ready-To-Use
Media


Irradiated
Media Bags

Pains

- Bottle handling
- Capacity/Volume
- Storage
- Autoclaving
- Weight of the bags
- Prewarming
- Manual preparation
- Risk of contamination
- Time consuming
- Multistep process




**Start the revolution with the
ReadyStream[®] system**



The ReadyStream® system

How does it work?




- Water



The ReadyStream® system


How does it work?



- Media



The ReadyStream® system




How does it work?

- Water




The ReadyStream® system

How does it work?



- Media
- Water





03 The benefits




Millipore®

Preparation, Separation,
Filtration & Monitoring Products

ReadyStream® Solutions

The benefits





**Prepare up to 100 L of media
in only about 20 min**

**(5 min of hands on and 15 minutes of
rehydration time)**



ReadyStream® Solutions

The benefits



10x concentrated media


Dispense up to 30 L with a 3 L bag

Dispense up to 100 L with a 10 L bag



ReadyStream® Solutions

The benefits





No more autoclaving



ReadyStream® Solutions

The benefits




Instant delivery of pre-heated media at different temperatures



ReadyStream® Solutions

The benefits




Dehydrated media bag

- Easy to carry and store
- Ready to be rehydrated



ReadyStream® Solutions

The benefits




5-days shelf life
for our media once
reconstituted



ReadyStream® Solutions

The benefits





ISO compliant with our GranuCult® culture media



ReadyStream® Solutions

The benefits




Small Footprint
saves space in your media kitchen




Solutions with the ReadyStream® system

The benefits




Flexible benchtop system that can also be used for a laminar flow hood





Improved sustainability







Q&A session



Millipore®

Preparation, Separation,
Filtration & Monitoring Products

**For more information or to book a
demo**

SigmaAldrich.com/ReadyStream-Demo



The vibrant M, Millipore, Milli-Q, SigmaAldrich, Supelco, ReadyStream and GranuCult are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2022 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.